


DERMALUX LED Face Mask

Built in clinics.
Perfected for home.





“
Every treatment
has taught us
something.
”

Welcome to *Dermalux*

As a British medical device manufacturer, Dermalux has pioneered professional LED Phototherapy since 2012, earning the trust of skincare professionals worldwide through clinical expertise, precision engineering and proven performance.

The Dermalux LED Face Mask brings this professional heritage into the home. Developed from the same clinically proven technology that powers thousands of professional treatments every day, it is tested to the same rigorous medical device standards as Dermalux professional systems. Bringing the benefits of professional LED Phototherapy into everyday skincare.



The Dermalux *Philosophy*

**Helping people feel more confident
through healthier skin.**

We believe confidence begins with healthier skin. Combining clinical expertise, science and precision engineering, we create technology trusted by skincare professionals worldwide to improve skin health and deliver safe, consistent and repeatable results.




Built on Biology. Engineered for *Results*.

Dermalux was founded on one simple belief: healthier skin comes from working with biology, not against it.

Powered by the science of photobiomodulation, Dermalux combines skin biology with precision engineering to optimise every treatment parameter, maximising the skin's biological response to deliver clinically proven results.

This biology-first philosophy defines every Dermalux device.



The Dermalux *Difference*

Professional LED Phototherapy isn't defined by a single specification. It's defined by the precision with which every treatment parameter works together.


Every Dermalux device begins with one simple question: **What does the skin need to respond?**

The answer defines *Dermalux DNA*. Our precision engineering philosophy.

From wavelength accuracy and treatment dose to treatment distance and light delivery, every treatment parameter is carefully engineered to work in synergy, optimising the skin's response to light for consistent, clinically proven results.

Not all LED devices are created equal.





Dermalux LED *Face Mask*

- Professional LED Phototherapy. Perfected for Home.
- Clinically Proven Skin Rejuvenation
- Dermalux Off-the-Skin Technology™
(Patent Pending)
- Three Clinically Proven Wavelengths
- One Expertly Optimised 10-Minute Treatment
- Dermalux Certified LEDs™
- Designed & Tested to Medical Device Standards





Dermalux *Off-The-Skin* Technology

More Coverage, *More Light Energy*.

Dermalux Off-the-Skin Technology™ positions the LEDs at the optimum treatment distance to maximise treatment coverage and deliver precision light energy for consistent, visible skin rejuvenation.

Up to 5x Greater Treatment Coverage*

Up to 2× More Light Energy*

**Based on independent testing against leading competitor brands*

Dermalux Off-The-Skin Technology (OTS) Patent Pending

Up to 5x Greater *Treatment Coverage**

Reaches more of the skin for a more complete treatment.



*96% Face Coverage**

Almost the entire treatment area.



Even Light Distribution



Extended Hairline Coverage


**Based on independent testing against leading competitor brands.*

**TYPICAL
ON-THE-SKIN MASKS**

Less than 20% face coverage*

10.6 j/cm²

at 0mm distance for 10 minutes




**DERMALUX
OFF-THE-SKIN TECHNOLOGY**

93% face coverage

21.2 j/cm²

at 10mm distance for 10 minutes



LOW COVERAGE



HIGH COVERAGE

**Tested against 6 competitor LED masks*

More coverage, *More light energy.*

Delivers more light where it matters most



Optimum Treatment Distance
Designed for efficient light delivery.



Reduced Hotspots
More even treatment.



Extended Hairline Coverage
Reliable treatment, every session.

**Based on independent testing against leading competitor brands.*


**TYPICAL
ON-THE-SKIN MASKS**



Hotspots



Coverage
Gaps




**DERMALUX
OFF-THE-SKIN TECHNOLOGY**



Even Light
Distribution



Optimised Light
Delivery



Three Clinically Proven *Wavelengths.*

One Complete Treatment.

RED 633nm | REJUVENATE


Smooths the appearance of fine lines and wrinkles
Improves skin tone and texture
Restores natural radiance

NEAR INFRARED 830nm | RESTORE

Supports skin repair and renewal
Reduces visible redness
Strengthens the skin barrier

DEEP NEAR INFRARED 1072nm | REVITALISE

Revives dull, tired-looking skin
Supports healthy microcirculation
Promotes skin vitality




One Mode *Expertly Optimised.*

Unlike multi-mode LED face masks, Dermalux delivers one expertly optimised 10-minute treatment where every critical treatment parameter has been engineered to work in synergy.

Clinically proven wavelengths. Optimised dose.
Optimised treatment time.

No modes. No settings. No guesswork. Just professional LED Phototherapy.





Designed & Tested Using *Medical Device Standards*

Built using the same rigorous quality standards as our professional systems.

Quality

Precision Engineered

Every treatment parameter optimised for maximum performance.

Safety

Medical Device Standards

Designed and tested to rigorous medical device standards.

Consistency

Consistent Performance

Reliable, repeatable treatment performance.

Confidence

Professional Heritage


Developed by a British medical device manufacturer.



Light That Fits *Your Life*

Comfortable. Convenient.
Consistent.

Professional LED Phototherapy, designed for everyday life



The Power of *Consistency*

Healthy skin is built over time.
Like regular exercise, the benefits
of LED Phototherapy are cumulative,
with each treatment building on the last.

Just 10 minutes a day.

Consistency turns clinically proven
light into lasting visible results.





Inside The *Box*

- Dermalux LED Face Mask
- Mask Controller
(connected to the mask)
- USB-C Charging Cable
- Protective Storage Pouch
- User Guide
- Warranty Card

Technical *Specification*

LED count:

144 LEDs (318 LED chips)

Optical output (irradiance):

36mw/cm² at 10mm

Treatment dose:

21.6 J/cm² at 10mm

Product dimensions:

190mm (W) x 190mm (D) x
270mm (H)

Controller dimensions:

40mm (W) x 45mm (D) x
130mm (H)

Cable length:

960mm

Battery recharge time:

2 hours 30 minutes

Treatment time:

10 minutes

Charger Type:

USB-C

Eye safety:

Integrated moulded eye
guards

Fitting type:

Adjustable head strap with
soft foam support

Controls:

Ergonomic one-button
controller

The Dermalux *Standard*

British by *Design*.
Professional by *Heritage*.
Trusted *Worldwide*.

Built on Biology. Backed by Science. Precision
engineered to deliver safe, consistent and
clinically proven results.





Discover the
Dermalux LED Face Mask




dermalux.com/products/dermalux-led-face-mask



dermaluxled



dermaluxled



Dermalux